



UBE2T Antibody, Biotin conjugated

Product Code	CSB-PA878841LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD8
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 T protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 T (EC 2.3.2.23) (Cell proliferation-inducing gene 50 protein) (E2 ubiquitin-conjugating enzyme T) (Ubiquitin carrier protein T) (Ubiquitin-protein ligase T), UBE2T
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBE2T
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.