



PIH1D3 Antibody, HRP conjugated

Product Code	CSB-PA878855LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQM4
Immunogen	Recombinant Human Protein PIH1D3 protein (1-214AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein PIH1D3 (PIH1 domain-containing protein 3) (Sarcoma antigen NY-SAR-97), PIH1D3, CXorf41
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PIH1D3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.