



CCDC87 Antibody

Product Code	CSB-PA878896LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVE4
Immunogen	Recombinant Human Coiled-coil domain-containing protein 87 protein (601-849AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 87, CCDC87
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCDC87