





SLC52A1 Antibody

Product Code	CSB-PA878904LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWF4
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 1 protein (212-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 1 (Porcine endogenous retrovirus A receptor 2) (PERV-A receptor 2) (Protein GPR172B) (Riboflavin transporter 1) (hRFT1), SLC52A1, GPR172B PAR2 RFT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC52A1


Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA878904LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA878904LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA878904LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.