



SLC52A1 Antibody, Biotin conjugated

Product Code	CSB-PA878904LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWF4
Immunogen	Recombinant Human Solute carrier family 52, riboflavin transporter, member 1 protein (212-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 52, riboflavin transporter, member 1 (Porcine endogenous retrovirus A receptor 2) (PERV-A receptor 2) (Protein GPR172B) (Riboflavin transporter 1) (hRFT1), SLC52A1, GPR172B PAR2 RFT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC52A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.