



GAR1 Antibody, Biotin conjugated

Product Code	CSB-PA878917LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY12
Immunogen	Recombinant Human H/ACA ribonucleoprotein complex subunit 1 protein (1-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	H/ACA ribonucleoprotein complex subunit 1 (Nucleolar protein family A member 1) (snRNP protein GAR1), GAR1, NOLA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.