



SMPD3 Antibody, Biotin conjugated

Product Code	CSB-PA878920LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NY59
Immunogen	Recombinant Human Sphingomyelin phosphodiesterase 3 protein (401-655AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sphingomyelin phosphodiesterase 3 (EC 3.1.4.12) (Neutral sphingomyelinase 2) (nSMase-2) (nSMase2) (Neutral sphingomyelinase II), SMPD3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SMPD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.