



MAP10 Antibody, FITC conjugated

Product Code	CSB-PA879069LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2G4
Immunogen	Recombinant Human Microtubule-associated protein 10 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein 10 (Microtubule regulator of 120 KDa), MAP10, KIAA1383 MTR120
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP10