



MAP10 Antibody, Biotin conjugated

Product Code	CSB-PA879069LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2G4
Immunogen	Recombinant Human Microtubule-associated protein 10 protein (401-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Microtubule-associated protein 10 (Microtubule regulator of 120 KDa), MAP10, KIAA1383 MTR120
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.