




Prok2 Antibody

Product Code	CSB-PA879446LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9QXU7
Immunogen	Recombinant Mouse Prokineticin-2 protein (27-128AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin-2 (PK2) (Protein Bv8 homolog), Prok2, Bv8
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Prok2

Image



Western blot

All lanes: Prok2 antibody at 3µg/ml + Mouse small intestine tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 15, 12, 9 kDa

Observed band size: 9 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.