



# Prok2 Antibody, Biotin conjugated

Product Code	CSB-PA879446LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9QXU7
Immunogen	Recombinant Mouse Prokineticin-2 protein (27-128AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prokineticin-2 (PK2) (Protein Bv8 homolog), Prok2, Bv8
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Prok2