



NIPSNAP3B Antibody, FITC conjugated

Product Code	CSB-PA880093LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS92
Immunogen	Recombinant Human Protein NipSnap homolog 3B protein (1-247AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NipSnap homolog 3B (NipSnap3B) (SNAP1), NIPSNAP3B, NIPSNAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NIPSNAP3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.