



NIPSNAP3B Antibody, Biotin conjugated

Product Code	CSB-PA880093LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS92
Immunogen	Recombinant Human Protein NipSnap homolog 3B protein (1-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein NipSnap homolog 3B (NipSnap3B) (SNAP1), NIPSNAP3B, NIPSNAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NIPSNAP3B