



TRIM55 Antibody, FITC conjugated

Product Code	CSB-PA880142LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYV6
Immunogen	Recombinant Human Tripartite motif-containing protein 55 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 55 (Muscle-specific RING finger protein 2) (MuRF-2) (MuRF2) (RING finger protein 29), TRIM55, MURF2 RNF29
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM55
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.