



Gp2 Antibody, HRP conjugated

Product Code	CSB-PA880511LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9D733
Immunogen	Recombinant Mouse Pancreatic secretory granule membrane major glycoprotein GP2 protein (22-531AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pancreatic secretory granule membrane major glycoprotein GP2 (Pancreatic zymogen granule membrane protein GP-2), Gp2
Immunogen Species	Mus musculus (Mouse)
Research Area	Microbiology
Target Names	Gp2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.