



SLC25A18 Antibody, FITC conjugated

Product Code	CSB-PA880945LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1K4
Immunogen	Recombinant Human Mitochondrial glutamate carrier 2 protein (246-287AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial glutamate carrier 2 (GC-2) (Glutamate/H(+) symporter 2) (Solute carrier family 25 member 18), SLC25A18, GC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.