



SLC25A18 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA880945LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H1K4 |
| Immunogen | Recombinant Human Mitochondrial glutamate carrier 2 protein (246-287AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Mitochondrial glutamate carrier 2 (GC-2) (Glutamate/H(+) symporter 2) (Solute carrier family 25 member 18), SLC25A18, GC2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | SLC25A18 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |