



SLC25A18 Antibody, Biotin conjugated

Product Code	CSB-PA880945LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1K4
Immunogen	Recombinant Human Mitochondrial glutamate carrier 2 protein (246-287AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial glutamate carrier 2 (GC-2) (Glutamate/H(+) symporter 2) (Solute carrier family 25 member 18), SLC25A18, GC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC25A18