



# SH3BGRL3 Antibody

<b>Product Code</b>	CSB-PA880952LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9H299
<b>Immunogen</b>	Recombinant Human SH3 domain-binding glutamic acid-rich-like protein 3 protein (2-93AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	SH3 domain-binding glutamic acid-rich-like protein 3 (SH3 domain-binding protein 1) (SH3BP-1), SH3BGRL3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	SH3BGRL3

## Image



### Western blot

All lanes: SH3BGRL3 antibody at 0.05μg/ml +


U87 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 11 kDa

Observed band size: 11 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA880952LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA880952LA01HU at dilution of 1:100