



SH3BGRL3 Antibody, Biotin conjugated

Product Code	CSB-PA880952LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H299
Immunogen	Recombinant Human SH3 domain-binding glutamic acid-rich-like protein 3 protein (2-93AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 domain-binding glutamic acid-rich-like protein 3 (SH3 domain-binding protein 1) (SH3BP-1), SH3BGRL3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SH3BGRL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.