



MRPL46 Antibody, HRP conjugated

Product Code	CSB-PA880959LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2W6
Immunogen	Recombinant Human 39S ribosomal protein L46, mitochondrial protein (1-279AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L46, mitochondrial (L46mt) (MRP-L46) (Mitochondrial large ribosomal subunit protein mL46) (P2ECSL), MRPL46, C15orf4 LIECG2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.