



GLIPR2 Antibody, Biotin conjugated

Product Code	CSB-PA880974LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4G4
Immunogen	Recombinant Human Golgi-associated plant pathogenesis-related protein 1 protein (2-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi-associated plant pathogenesis-related protein 1 (GAPR-1) (Golgi-associated PR-1 protein) (Glioma pathogenesis-related protein 2) (GliPR 2), GLIPR2, C9orf19 GAPR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GLIPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.