



# GLIPR2 Antibody, Biotin conjugated

Product Code	CSB-PA880974LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4G4
Immunogen	Recombinant Human Golgi-associated plant pathogenesis-related protein 1 protein (2-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi-associated plant pathogenesis-related protein 1 (GAPR-1) (Golgi-associated PR-1 protein) (Glioma pathogenesis-related protein 2) (GliPR 2), GLIPR2, C9orf19 GAPR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GLIPR2