



# RABEP2 Antibody, FITC conjugated

Product Code	CSB-PA880980LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5N1
Immunogen	Recombinant Human Rab GTPase-binding effector protein 2 protein (301-533AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab GTPase-binding effector protein 2 (Rabaptin-5beta), RABEP2, RABPT5B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RABEP2