



C17orf75 Antibody

Product Code	CSB-PA881009LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAS0
Immunogen	Recombinant Human Protein Njmu-R1 protein (201-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Njmu-R1, C17orf75
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C17orf75
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.