



SPC25 Antibody, Biotin conjugated

Product Code	CSB-PA881019LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBM1
Immunogen	Recombinant Human Kinetochores protein Spc25 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetochores protein Spc25 (hSpc25), SPC25, SPBC25
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SPC25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.