




VPS45 Antibody

Product Code	CSB-PA882092LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRW7
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 45 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 45 (h-VPS45) (hVps45), VPS45, VPS45A VPS45B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS45

Image



Western blot

All lanes: VPS45 antibody at 3μg/ml

Lane 1: Mouse brain tissue


Lane 2: U87 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 66, 62 kDa

Observed band size: 41 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA882092LA01HU at dilution of 1:100