



VPS45 Antibody, HRP conjugated

Product Code	CSB-PA882092LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRW7
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 45 protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 45 (h-VPS45) (hVps45), VPS45, VPS45A VPS45B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS45
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.