



ECHDC1 Antibody

Product Code	CSB-PA882103LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTX5
Immunogen	Recombinant Human Ethylmalonyl-CoA decarboxylase protein (1-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ethylmalonyl-CoA decarboxylase (EC 4.1.1.94) (Enoyl-CoA hydratase domain-containing protein 1) (Methylmalonyl-CoA decarboxylase) (MMCD) (EC 4.1.1.41), ECHDC1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ECHDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.