



GID8 Antibody, HRP conjugated

Product Code	CSB-PA882124LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWU2
Immunogen	Recombinant Human Glucose-induced degradation protein 8 homolog protein (1-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-induced degradation protein 8 homolog (Two hybrid-associated protein 1 with RanBPM) (Twa1), GID8, C20orf11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GID8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.