



LYAR Antibody, Biotin conjugated

Product Code	CSB-PA882130LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX58
Immunogen	Recombinant Human Cell growth-regulating nucleolar protein (1-379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell growth-regulating nucleolar protein, LYAR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LYAR