



ELOVL5 Antibody, Biotin conjugated

Product Code	CSB-PA882144LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYP7
Immunogen	Recombinant Human Elongation of very long chain fatty acids protein 5 protein (247-299AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Elongation of very long chain fatty acids protein 5 (EC 2.3.1.199) (3-keto acyl-CoA synthase ELOVL5) (ELOVL fatty acid elongase 5) (ELOVL FA elongase 5) (Fatty acid elongase 1) (hELO1) (Very long chain 3-ketoacyl-CoA synthase 5) (Very long chain 3-oxoacyl-CoA synthase 5), ELOVL5, ELOVL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ELOVL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.