



PKD2L2 Antibody, HRP conjugated

Product Code	CSB-PA882158LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZM6
Immunogen	Recombinant Human Polycystic kidney disease 2-like 2 protein (390-523AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycystic kidney disease 2-like 2 protein (Polycystin-2L2) (Polycystin-L2), PKD2L2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PKD2L2