



CCDC59 Antibody, FITC conjugated

Product Code	CSB-PA882165LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P031
Immunogen	Recombinant Human Thyroid transcription factor 1-associated protein 26 protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroid transcription factor 1-associated protein 26 (TTF-1-associated protein 26) (Coiled-coil domain-containing protein 59) (TTF-1-associated protein BR2), CCDC59, BR22 TAP26
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	CCDC59
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.