



SEC14L4 Antibody, Biotin conjugated

Product Code	CSB-PA883374LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UDX3
Immunogen	Recombinant Human SEC14-like protein 4 protein (201-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SEC14-like protein 4 (Tocopherol-associated protein 3), SEC14L4, TAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC14L4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.