



FBXO40 Antibody, HRP conjugated

Product Code	CSB-PA883385LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UH90
Immunogen	Recombinant Human F-box only protein 40 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 40 (Muscle disease-related protein), FBXO40, FBX40 KIAA1195
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO40
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.