



CLEC2D Antibody, FITC conjugated

Product Code	CSB-PA883401LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHP7
Immunogen	Recombinant Human C-type lectin domain family 2 member D protein (60-154AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 2 member D (Lectin-like NK cell receptor) (Lectin-like transcript 1) (LLT-1) (Osteoclast inhibitory lectin), CLEC2D, CLAX LLT1 OCIL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CLEC2D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.