



CLEC2D Antibody, Biotin conjugated

Product Code	CSB-PA883401LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHP7
Immunogen	Recombinant Human C-type lectin domain family 2 member D protein (60-154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 2 member D (Lectin-like NK cell receptor) (Lectin-like transcript 1) (LLT-1) (Osteoclast inhibitory lectin), CLEC2D, CLAX LLT1 OCIL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CLEC2D