



# MRM2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA883420LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UI43
<b>Immunogen</b>	Recombinant Human rRNA methyltransferase 2, mitochondrial protein (19-246AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	rRNA methyltransferase 2, mitochondrial (EC 2.1.1.-) (16S rRNA (uridine(1369)-2'-O)-methyltransferase) (16S rRNA [Um1369] 2'-O-methyltransferase) (Protein ftsJ homolog 2), MRM2, FJH1 FTSJ2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	MRM2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.