



MRM2 Antibody, Biotin conjugated

Product Code	CSB-PA883420LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI43
Immunogen	Recombinant Human rRNA methyltransferase 2, mitochondrial protein (19-246AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	rRNA methyltransferase 2, mitochondrial (EC 2.1.1.-) (16S rRNA (uridine(1369)-2'-O)-methyltransferase) (16S rRNA [Um1369] 2'-O-methyltransferase) (Protein ftsJ homolog 2), MRM2, FJH1 FTSJ2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRM2