



ST6GALNAC2 Antibody, Biotin conjugated

Product Code	CSB-PA883452LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJ37
Immunogen	Recombinant Human Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 protein (101-374AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-N-acetylgalactosaminide alpha-2,6-sialyltransferase 2 (EC 2.4.99.-) (GalNAc alpha-2,6-sialyltransferase II) (ST6GalNAc II) (ST6GalNAcII) (SThM) (Sialyltransferase 7B) (SIAT7-B), ST6GALNAC2, SIAT7B SIATL1 STHM
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ST6GALNAC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.