



MS4A8 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA883629LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BY19 |
| Immunogen | Recombinant Human Membrane-spanning 4-domains subfamily A member 8 protein (20-74AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Membrane-spanning 4-domains subfamily A member 8 (Four-span transmembrane protein 4) (Membrane-spanning 4-domains subfamily A member 8B), MS4A8, 4SPAN4 MS4A8B |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | MS4A8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |