



GSDMC Antibody, FITC conjugated

Product Code	CSB-PA883633LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYG8
Immunogen	Recombinant Human Gasdermin-C protein (1-508AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gasdermin-C (Melanoma-derived leucine zipper-containing extranuclear factor), GSDMC, MLZE
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GSDMC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.