



OSBPL6 Antibody, FITC conjugated

Product Code	CSB-PA883645LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZF3
Immunogen	Recombinant Human Oxysterol-binding protein-related protein 6 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterol-binding protein-related protein 6 (ORP-6) (OSBP-related protein 6), OSBPL6, ORP6
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OSBPL6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.