




Snrnp35 Antibody

Product Code	CSB-PA883986LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9D384
Immunogen	Recombinant Mouse U11/U12 small nuclear ribonucleoprotein 35 kDa protein (1-244AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U11/U12 small nuclear ribonucleoprotein 35 kDa protein (U11/U12 snRNP 35 kDa protein) (U1 snRNP-binding protein homolog), Snrnp35, U1snrnpbp
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Snrnp35

Image



Western blot

All lanes: Snrnp35 antibody at 12µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 29 kDa
 Observed band size: 29 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.