



RIBC2 Antibody, Biotin conjugated

Product Code	CSB-PA884469LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4K1
Immunogen	Recombinant Human RIB43A-like with coiled-coils protein 2 protein (1-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RIB43A-like with coiled-coils protein 2, RIBC2, C22orf11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RIBC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.