





VSIR Antibody

Product Code	CSB-PA884484LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7M9
Immunogen	Recombinant Human V-type immunoglobulin domain-containing suppressor of T-cell activation protein (37-194AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	V-type immunoglobulin domain-containing suppressor of T-cell activation (Platelet receptor Gi24) (Stress-induced secreted protein-1) (Sisp-1) (V-set domain-containing immunoregulatory receptor) (V-set immunoregulatory receptor), VSIR, C10orf54 SISP1 VISTA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VSIR

Image



Western blot
All lanes: C10orf54 antibody at 4μg/ml + Mouse liver tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 34 kDa
Observed band size: 34 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA884484LA01HU at dilution of 1:100