



VSIR Antibody, HRP conjugated

Product Code	CSB-PA884484LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7M9
Immunogen	Recombinant Human V-type immunoglobulin domain-containing suppressor of T-cell activation protein (37-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	V-type immunoglobulin domain-containing suppressor of T-cell activation (Platelet receptor Gi24) (Stress-induced secreted protein-1) (Sisp-1) (V-set domain-containing immunoregulatory receptor) (V-set immunoregulatory receptor), VSIR, C10orf54 SISP1 VISTA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	VSIR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.