



NMUR1 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA884509LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9HB89 |
| Immunogen | Recombinant Human Neuromedin-U receptor 1 protein (360-426AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Neuromedin-U receptor 1 (NMU-R1) (G-protein coupled receptor 66) (G-protein coupled receptor FM-3), NMUR1, GPR66 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | NMUR1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |