



NMUR1 Antibody, Biotin conjugated

Product Code	CSB-PA884509LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HB89
Immunogen	Recombinant Human Neuromedin-U receptor 1 protein (360-426AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuromedin-U receptor 1 (NMU-R1) (G-protein coupled receptor 66) (G-protein coupled receptor FM-3), NMUR1, GPR66
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NMUR1